

Welcome

While at the collection site, please practice social distancing and wear a mask except when providing saliva. The Shield T3 test is non-invasive and shouldn't take more than 3 minutes. Follow the guidelines below and any instructions you receive from site staff.



ShieldT3
Health
TARGET • TEST • TELL

Prep



Do not eat, drink, smoke, or brush your teeth within an hour of your test.



Sanitize your hands.



Retrieve a vial and funnel.

Provide saliva




Remove the cap and hold on to it.



Attach the funnel securely to the vial.



Deposit 1 to 1.5 ml of saliva into the funnel.



Remove the funnel and discard it in a medical waste container.



Screw the cap on the vial and tighten so the vial is completely sealed.

Scan



Sanitize the vial and your hands with a sanitizing wipe.



Write your birth month and day on the vial.



Retrieve a bag. Put your vial in the bag, and seal.



Present the vial and QR code for scanning.



Place the scanned sample where directed.